

Anti-Kir2DL4 Polyclonal Antibody For Wb, Elisa - Q99706

Item Number: CM0012909



Introduction

KIR2DL4 Rabbit Polyclonal Antibody (Q99706) validated for WB, ELISA. Detects human, mouse, rat KIR2DL4, a receptor for non-classical HLA-G. Ideal for studying NK cell signaling and maternal-fetal immunology.

[Learn More](#)

Parameter	Value
Product Name	KIR2DL4 Rabbit pAb
Remarks/Alias	G9P; CD158D; KIR103; KIR-2DL4; KIR103AS; KIR-103AS; KIR2DL4
Species	Human
Gene ID (Human)	3805
Gene ID	3805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-242 of human KIR2DL4 (NP_002246.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99706
Protein Weight	41kDa
Shipping	Ice bag