

Anti-Klf10 Polyclonal Antibody For Wb, Elisa - Q13118

Item Number: CM0018610



Introduction

Rabbit polyclonal antibody targeting human KLF10 (TIEG1). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for studies on transcriptional regulation, circadian rhythms, and metabolism.

[Learn More](#)

Parameter	Value
Product Name	KLF10 Rabbit pAb
Remarks/Alias	EGRA; TIEG; TIEG1; EGR-alpha; KLF10
Species	Human
Gene ID (Human)	7071
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13118
Protein Weight	53kDa
Shipping	Ice bag