

Anti-Klf8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95600

Item Number: CM0023084

Introduction

Recombinant KLF8 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for studying transcription regulation, cell cycle progression, and KLF8-related pathways.

[Learn More](#)

Parameter	Value
Product Name	KLF8 Rabbit pAb
Remarks/Alias	BKLF3; ZNF741; KLF8
Species	Human
Gene ID (Human)	11279
Gene ID	11279
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-280 of human KLF8 (NP_009181.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95600
Protein Weight	39kDa
Shipping	Ice bag