


Anti-Klkb1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P03952

Item Number: CM0029526



Introduction

High-quality rabbit polyclonal antibody against human KLKB1 (plasma kallikrein). Validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for coagulation, bradykinin, and renin-angiotensin research.

[Learn More](#)

Parameter	Value
Product Name	KLKB1 Rabbit pAb
Remarks/Alias	PKK; PPK; KLK3; PKKD; KLKB1
Species	Human
Gene ID (Human)	3818
Gene ID	3818
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human KLKB1 (NP_000883.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03952
Protein Weight	71kDa
Shipping	Ice bag