


Anti-Ku80 Monoclonal Antibody For Elisa, If/Icc, Ihc-P, Ip, Wb - P13010

Item Number: CM0005433



Introduction

Recombinant rabbit monoclonal antibody against human Ku80 (XRCC5). Validated for WB, IP, IF/ICC, IHC-P, and ELISA. Suitable for NHEJ, DNA repair, and telomere maintenance research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Ku80 Rabbit mAb
Remarks/Alias	KU80; KUB2; Ku86; NFIV; KARP1; KARP-1
Species	Human
Gene ID (Human)	7520
Gene ID	7520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 543-732 of human Ku80 (NP_066964.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P13010
Protein Weight	83 kDa
Shipping	Ice bag