


Anti-Lactb Polyclonal Antibody For Wb And Elisa - P83111

Item Number: CM0030165



Introduction

Rabbit polyclonal antibody targeting human LACTB, a mitochondrial protease and tumor suppressor. Recommended for WB and ELISA; cross-reacts with human and mouse. Ideal for mitochondrial biology and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	LACTB Rabbit pAb
Remarks/Alias	G24; MRPL56; LACTB
Species	Human
GeneID (Human)	114294
GeneID	114294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-290 of human LACTB (NP_116246.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P83111
Protein Weight	61kDa
Shipping	Ice bag