

Anti-Lair2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q6Iss4

Item Number: CM0014654

Introduction

High-quality rabbit polyclonal antibody against human LAIR2 (CD306), suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for leukocyte and immune receptor studies.

[Learn More](#)

Parameter	Value
Product Name	LAIR2 Rabbit pAb
Remarks/Alias	CD306; LAIR2
Species	Human
GeneID (Human)	3904
GeneID	3904
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-120 of human LAIR2 (NP_002279.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ISS4
Protein Weight	16kDa
Shipping	Ice bag