

Anti-Lamin A/C Polyclonal Antibody For Wb, If/Icc, Ihc-P, Chip, Elisa - P02545

Item Number: CM0016934



Introduction

Rabbit polyclonal antibody targeting human Lamin A/C, suitable for Western blot, immunocytochemistry, immunohistochemistry, chromatin immunoprecipitation, and ELISA. Recognizes human, mouse, and rat species. Ideal for studies of nuclear lamina structure, chromosome dynamics, and laminopathy mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Lamin A/C Rabbit pAb
Remarks/Alias	FPL; IDC; LFP; CDDC; EMD2; FPLD; HGPS; LDP1; LMN1; LMNC; MADA; PRO1; CDCD1; CMD1A; FPLD2; LMNL1; CMT2B1; LGMD1B; Lamin A/C
Species	Human
Gene ID (Human)	4000
Gene ID	4000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 403-572 of human Lamin A/C (NP_733821.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, Chip, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02545
Protein Weight	74 kDa
Shipping	Ice bag