

Anti-Langerin/Cd207 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uj71

Item Number: CM0026597



Introduction

Rabbit polyclonal antibody against Langerin/CD207. Validated for Western blotting and ELISA; reacts with human, mouse, and rat. Ideal for studying Langerin function in immune cells and antigen presentation.

[Learn More](#)

Parameter	Value
Product Name	Langerin/CD207 Rabbit pAb
Remarks/Alias	CLEC4K; Langerin/CD207
Species	Human
GeneID (Human)	50489
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-200 of human Langerin/CD207 (NP_056532.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJ71
Protein Weight	37kDa
Shipping	Ice bag