

Anti-Larp1B Polyclonal Antibody For Wb, Ihc-P, Elisa - Q659C4

Item Number: CM0016982

Introduction

Rabbit polyclonal antibody against human LARP1B. Validated for Western blot, IHC-P, and ELISA. Detects LARP1B (~105 kDa) in human, mouse, and rat samples. Ideal for RNA-binding protein research and assay development.

[Learn More](#)

Parameter	Value
Product Name	LARP1B Rabbit pAb
Remarks/Alias	LARP2; LARP1B
Species	Human
GeneID (Human)	55132
GeneID	55132
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-320 of human LARP1B (NP_060548.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q659C4
Protein Weight	105kDa
Shipping	Ice bag