

Anti-Lefty1/2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - 075610 / 000292

Item Number: CM0004170

Introduction

Rabbit monoclonal antibody against human LEFTY1/2 proteins, validated for WB, IHC-P, and ELISA. Ideal for detecting left-right determination factors in developmental research across human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	LEFTY1/2 Rabbit mAb
Remarks/Alias	LEFTB; LEFTYB; LEFTY1/2
Species	Human
Gene ID (Human)	10637/7044
Gene ID	10637/7044
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 267-366 of human LEFTY1/2 (075610).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075610/000292
Protein Weight	34.44kDa
Shipping	Ice bag