

## Anti-Lepre1 (Gros1, Oi8) Polyclonal Antibody For Wb, Elisa - Q32P28

Item Number: CM0014754



### Introduction

Rabbit polyclonal antibody targeting human Prolyl 3-hydroxylase 1 (LEPRE1/GROS1/OI8), suitable for Western blot and ELISA applications. Useful for collagen biosynthesis and skeletal disease research.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | LEPRE1 Rabbit pAb (LEPRE1)                                                                                           |
| Remarks/Alias    | OI8; GROS1; LEPRE1                                                                                                   |
| Species          | Human                                                                                                                |
| Gene ID          | 64175                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 497-736 of human LEPRE1 (NP_071751.3). |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human                                                                                                                |
| SWISS            | Q32P28                                                                                                               |
| Protein Weight   | 83kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |