

Anti-Leptin Receptor Monoclonal Antibody For Wb, Ihc-P, Elisa - P48357

Item Number: CM0003721



Introduction

Leptin Receptor Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Cross-reacts with Human, Mouse, Rat. Detects ~132 kDa LEPR, a key regulator of appetite and energy balance.

[Learn More](#)

Parameter	Value
Product Name	Leptin Receptor Rabbit mAb
Remarks/Alias	OBR; OB-R; CD295; LEP-R; LEPRD; Leptin Receptor
Species	Human
GeneID (Human)	3953
GeneID	3953
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1000-1100 of human Leptin Receptor (P48357).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48357
Protein Weight	132kDa
Shipping	Ice bag