

Anti-Lf1 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Taq6

Item Number: CM0020562

Introduction

Rabbit polyclonal antibody specific to *Oryza sativa* LF1, a homeobox-leucine zipper transcription factor (HOX10/OsHB1). Validated for WB and ELISA, suitable for gene regulation and plant development research.

[Learn More](#)

Parameter	Value
Product Name	LF1 Rabbit pAb
Remarks/Alias	hb1; HOX10; OsHB1; OsHox10; Osj_008792; LF1
Species	<i>Oryza sativa</i>
GeneID	4331345
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 400-500 of <i>Oryza sativa</i> LF1. (Q6TAQ6).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q6TAQ6
Protein Weight	92kDa
Shipping	Ice bag