

Anti-Lgi1 Rabbit Monoclonal Antibody For Wb And Elisa - O95970

Item Number: CM0003521



Introduction

Anti-LGI1 rabbit monoclonal antibody for Western blot (WB) and ELISA, targeting human LGI1 (O95970) with cross-reactivity to mouse and rat. Suitable for neuroscience research on synaptic transmission and epilepsy mechanisms.

[Learn More](#)

Parameter	Value
Product Name	LGI1 Rabbit mAb
Remarks/Alias	EPT; ETL1; ADLTE; ADPAEF; ADPEAF; IB1099; EPITEMPIN; LGI1
Species	Human
GeneID (Human)	9211
GeneID	9211
Immunogen	Synthetic peptide (A synthetic peptide corresponding to a sequence within amino acids 150-250 of human LGI1 (O95970))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O95970
Protein Weight	64kDa
Shipping	Ice bag