

Anti-Lgr4 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Bxb1

Item Number: CM0014533



Introduction

LGR4 Rabbit Polyclonal Antibody validated for WB, IF-P, ELISA. Reacts with human, mouse, rat LGR4, a Wnt signaling receptor crucial for organ development, stem cell maintenance, and bone remodeling. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	LGR4 Rabbit pAb
Remarks/Alias	DPSL; GPR48; BNMD17; LGR4
Species	Human
Gene ID (Human)	55366
Gene ID	55366
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-540 of human LGR4 (NP_060960.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXB1
Protein Weight	104kDa
Shipping	Ice bag