

## Anti-Lhx1 Polyclonal Antibody For Wb, Elisa - P48742

Item Number: CM0021512

### Introduction

Rabbit polyclonal antibody targeting human LHX1 (LIM1), validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Suitable for developmental biology and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | LHX1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | LIM1; LIM-1; LHX1                                                                                                 |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 3975                                                                                                              |
| GeneID           | 3975                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 120-180 of human LHX1 (NP_005559.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | P48742                                                                                                            |
| Protein Weight   | 45kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |