

Anti-Lif Rabbit Polyclonal Antibody For Wb, Elisa - P09056

Item Number: CM0019993



Introduction

CamelBio offers a validated rabbit polyclonal antibody against mouse leukemia inhibitory factor (LIF) for WB and ELISA, ideal for studying embryonic stem cell self-renewal and implantation.

[Learn More](#)

Parameter	Value
Product Name	Lif Rabbit pAb
Species	Mouse
GeneID (Human)	3976
GeneID	16878
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-203 of mouse Lif (NP_032527.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P09056
Protein Weight	22kDa
Shipping	Ice bag