


Anti-Lipa Polyclonal Antibody For Wb, Elisa - P38571

Item Number: CM0030089



Introduction

Rabbit polyclonal antibody against human LIPA (lysosomal acid lipase), validated for WB and ELISA. Cross-reacts with rat. Suitable for research in lipid metabolism and lysosomal storage disorders.

[Learn More](#)

Parameter	Value
Product Name	LIPA Rabbit pAb
Remarks/Alias	LAL; CESD; LIPA
Species	Human
GeneID (Human)	3988
GeneID	3988
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-399 of human LIPA (NP_001121077.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P38571
Protein Weight	45kDa
Shipping	Ice bag