

Anti-Liph Rabbit Polyclonal Antibody For Wb & Elisa - Q8Wwy8

Item Number: CM0018748

Introduction

Rabbit polyclonal antibody targeting human LIPH (lipase member H), a phosphatidic acid-selective phospholipase A1. Suitable for Western blotting (WB) and ELISA; cross-reacts with human, mouse, and rat. Supports lipid metabolism and hair biology studies.

[Learn More](#)

Parameter	Value
Product Name	LIPH Rabbit pAb
Remarks/Alias	AH; LAH2; ARWH2; HYPT7; LPDLR; PLA1B; mPA-PLA1; LIPH
Species	Human
GeneID (Human)	200879
GeneID	200879
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human LIPH (NP_640341.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WWY8
Protein Weight	51kDa
Shipping	Ice bag