


Anti-Lipocalin-2/NGal Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P80188 / P11672

Item Number: CM0011812



Introduction

Rabbit polymorphic monoclonal antibody against human/mouse Lipocalin-2/NGAL. Validated for WB, IF-P, IHC-P, and ELISA. Suitable for apoptosis, innate immunity, and iron metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Lipocalin-2/NGAL Rabbit PolymAb®
Remarks/Alias	p25; 24p3; MSFI; NGAL; Lipocalin-2/NGAL; AW212229; NRL; Sip24
Species	Human/Mouse
GeneID (Human)	3934
GeneID	3934/16819
Immunogen	Recombinant protein of human/mouse Lipocalin-2/NGAL
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P80188/P11672
Protein Weight	23kDa
Shipping	Ice bag