

Anti-Lpcat4 Polyclonal Antibody For Wb, Elisa - Q643R3

Item Number: CM0018277

Introduction

Rabbit polyclonal antibody against human LPCAT4, reactive with human and mouse. Ideal for WB and ELISA applications. LPCAT4 is a lysophospholipid acyltransferase involved in phospholipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	LPCAT4 Rabbit pAb
Remarks/Alias	AYTL3; AGPAT7; LPEAT2; LPLAT10; LPAAT-eta; LPCAT4
Species	Human
GeneID (Human)	254531
GeneID	254531
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 345-524 of human LPCAT4 (NP_705841.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q643R3
Protein Weight	57kDa
Shipping	Ice bag