


Anti-Lsm3 Rabbit Polyclonal Antibody For Wb, Elisa - P62310

Item Number: CM0016192



Introduction

LSM3 rabbit polyclonal antibody validated for Western blot (WB) and ELISA applications, detecting human and mouse U6 snRNA-associated Sm-like protein LSM3, a spliceosome component.

[Learn More](#)

Parameter	Value
Product Name	LSM3 Rabbit pAb
Remarks/Alias	SMX4; USS2; YLR438C; LSM3
Species	Human
Gene ID (Human)	27258
Gene ID	27258
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-102 of human LSM3 (NP_055278.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P62310
Protein Weight	12kDa
Shipping	Ice bag