


Anti-Lss Rabbit Polyclonal Antibody For Wb And Elisa - P48449

Item Number: CM0025704



Introduction

Rabbit polyclonal antibody targeting human LSS, the key enzyme in cholesterol biosynthesis, validated for WB and ELISA, cross-reacts with mouse. Ideal for research on lanosterol synthesis.

[Learn More](#)

Parameter	Value
Product Name	LSS Rabbit pAb
Remarks/Alias	OSC, APMR4, HYPT14, CTRCT44, LSS
Species	Human
Gene ID (Human)	4047
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human LSS (NP_002331.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48449
Protein Weight	83kDa
Shipping	Ice bag