

Anti-Lxr Beta/Ner Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55055

Item Number: CM0023103



Introduction

LXR beta/NER rabbit polyclonal antibody, validated for WB, IF/ICC, ELISA, reacts with human, mouse, rat. Targets oxysterols receptor LXR-beta for nuclear receptor and cholesterol metabolism research.

[Learn More](#)

Parameter	Value
Product Name	LXR beta/NER Rabbit pAb
Remarks/Alias	NER; UNR; LXRβ; LXR-b; NER-I; RIP15; NR1H2
Species	Human
GeneID (Human)	7376
GeneID	7376
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 180-260 of human LXR beta/NER (NP_009052.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55055
Protein Weight	51kDa
Shipping	Ice bag