

Anti-Ly86 Rabbit Polyclonal Antibody For Wb, Elisa - O95711

Item Number: CM0019856



Introduction

High-quality rabbit polyclonal antibody targeting human LY86 (MD1). Validated for WB and ELISA, with cross-reactivity to mouse. Ideal for innate immunity and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	LY86 Rabbit pAb
Remarks/Alias	MD1; MD-1; MMD-1; dJ80N2.1; LY86
Species	Human
GeneID (Human)	9450
GeneID	9450
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95711
Protein Weight	18kDa
Shipping	Ice bag