


Anti-Lysrs/Kars Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q15046

Item Number: CM0009758



Introduction

Anti-LysRS/KARS Rabbit Monoclonal Antibody for detection of Lysyl-tRNA synthetase in human, mouse, rat samples using WB, IHC-P, IF/ICC, IP, ELISA. Ideal for IVD raw material and research applications.

[Learn More](#)

Parameter	Value
Product Name	LysRS/KARS Rabbit mAb
Remarks/Alias	KRS; KARS; KARS2; LEPID; CMTRIB; DEAPLE; DFN89; LysRS/KARS
Species	Human
GeneID (Human)	3735
GeneID	3735
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15046
Protein Weight	68kDa
Shipping	Ice bag