


Anti-Lyve1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y5Y7

Item Number: CM0025601



Introduction

Rabbit polyclonal antibody recognizing LYVE1 (lymphatic vessel endothelial hyaluronic acid receptor 1) for Western blot and ELISA applications. Reacts with human and mouse samples, ideal for lymphatic vessel and hyaluronan research.

[Learn More](#)

Parameter	Specification
Product Name	LYVE1 Rabbit pAb
Remarks/Alias	HAR; XLKD1; LYVE-1; CRSBP-1; LYVE1
Species	Human
GeneID (Human)	10894
GeneID	10894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-238 of human LYVE1 (NP_006682.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y5Y7
Protein Weight	35kDa
Shipping	Ice bag