


Anti-MacroH2A.1 Monoclonal Antibody For Wb, IHC-P, IF/ICC, Elisa - O75367

Item Number: CM0008132



Introduction

Recombinant rabbit monoclonal macroH2A.1 antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, and rat macroH2A.1, a variant histone involved in transcriptional repression and NAD⁺ metabolism. Ideal for chromatin research.

[Learn More](#)

Parameter	Value
Product Name	macroH2A.1 Rabbit mAb
Remarks/Alias	H2A.y; H2A/y; H2AFY; mH2A1; H2AF12M; MACROH2A1.1; macroH2A1.2; macroH2A.1
Species	Human
GeneID (human)	9555
GeneID	9555
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human macroH2A.1 (O75367)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75367
Protein Weight	39kDa
Shipping	Ice bag