

Anti-Mad2/Mad2L1 Rabbit Monoclonal Antibody For Wb, Elisa - Q13257

Item Number: CM0003058



Introduction

Anti-MAD2/MAD2L1 rabbit monoclonal antibody for WB and ELISA applications. Designed to detect human MAD2A, a critical spindle assembly checkpoint protein. Suitable for cell cycle and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	MAD2/MAD2L1 Rabbit mAb
Remarks/Alias	MAD2; HSMAD2; MAD2/MAD2L1
Species	Human
Gene ID (Human)	4085
Gene ID	4085
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 70-150 of human MAD2/MAD2L1 (Q13257).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13257
Protein Weight	24kDa
Shipping	Ice bag