

Anti-Madcam-1 Polyclonal Antibody For Wb, Elisa - Q61826

Item Number: CM0020470

Introduction

Rabbit polyclonal antibody against mouse MAdCAM-1, validated for Western blot and ELISA. Targets mucosal addressin cell adhesion molecule 1, a key regulator of lymphocyte homing to mucosal tissues.

[Learn More](#)

Parameter	Value
Product Name	MAdCAM-1 Rabbit pAb
Remarks/Alias	MAdCAM-1
Species	Mouse
GeneID	17123
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q61826
Protein Weight	43kDa
Shipping	Ice bag