

Anti-Man1C1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nr34

Item Number: CM0015710

Introduction

Rabbit polyclonal antibody targeting human MAN1C1, a Golgi alpha-1,2-mannosidase involved in N-glycan processing. Validated for western blot (WB) and ELISA, with cross-reactivity to mouse. Ideal for glycobiology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MAN1C1 Rabbit pAb
Remarks/Alias	HMIC; MAN1C; MAN1A3; pp6318; MAN1C1
Species	Human
GeneID (Human)	57134
GeneID	57134
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 44-240 of human MAN1C1 (NP_065112.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NR34
Protein Weight	71kDa
Shipping	Ice bag