


Anti-Map1Lc3B Polyclonal Antibody For Wb, Elisa - Q9Gzq8

Item Number: CM0020949



Introduction

MAP1LC3B Rabbit Polyclonal Antibody validated for WB and ELISA with human and rat cross-reactivity. Essential for studying autophagy, mitophagy, and autophagosome formation in IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	MAP1LC3B Rabbit pAb
Remarks/Alias	LC3B; ATG8F; MAP1LC3B-a; MAP1A/1BLC3; 3B
Species	Human
Gene ID (Human)	81631
Immunogen	Synthetic peptide within amino acids 1-50 of human MAP1LC3B (NP_073729.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9GZQ8
Protein Weight	14kDa
Shipping	Ice bag