


Anti-Mars Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P56192

Item Number: CM0027328



Introduction

High-quality rabbit polyclonal antibody against human MARS (methionine-tRNA ligase), validated for WB, IHC-P, and ELISA. Detects human and rat MARS, ideal for protein synthesis and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	MARS Rabbit pAb
Remarks/Alias	MRS; ILLD; MARS; TTD9; CMT2U; ILFS2; METRS; MTRNS; SPG70
Species	Human
Gene ID (Human)	4141
Gene ID	4141
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human MARS (NP_004981.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P56192
Protein Weight	101kDa
Shipping	Ice bag