


Anti-Marveld2 (Tricellulin) Rabbit Polyclonal Antibody For Wb, Elisa - Q8N4S9

Item Number: CM0012924



Introduction

Rabbit polyclonal antibody targeting human MARVELD2 (tricellulin), validated for WB and ELISA. Reacts with human, mouse, and rat samples. Ideal for tight junction and deafness research.

[Learn More](#)

Parameter	Value
Product Name	MARVELD2 Rabbit pAb
Remarks/Alias	Tric; DFN49; MARVD2; MRVLC2; MARVELD2
Species	Human
GeneID (Human)	153562
GeneID	153562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human MARVELD2 (NP_001033692.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N4S9
Protein Weight	64kDa
Shipping	Ice bag