


Anti-Mc2R Polyclonal Antibody For Wb And Elisa - Q01718

Item Number: CM0029901



Introduction

High-quality MC2R rabbit polyclonal antibody validated for WB and ELISA. Detects endogenous MC2R in human and mouse samples. Ideal for investigating adrenocortical function, cortisol release, and GPCR signaling.

[Learn More](#)

Parameter	Value
Product Name	MC2R Rabbit pAb
Remarks/Alias	ACTHR; MC2R
Species	Human
GeneID (Human)	4158
GeneID	4158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human MC2R (NP_000520.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q01718
Protein Weight	34kDa
Shipping	Ice bag