

Anti-Mchr1 Polyclonal Antibody For Wb And Elisa - Q99705

Item Number: CM0021737



Introduction

Rabbit anti-MCHR1 polyclonal antibody for WB and ELISA. Cross-reacts with human, mouse, rat. MCHR1 is a GPCR regulating feeding, energy balance, and sleep. Ideal for neuroscience and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	MCHR1 Rabbit pAb
Remarks/Alias	SLC1; GPR24; MCH1R; SLC-1; MCH-1R; MCHR1
Species	Human
GeneID (Human)	2847
GeneID	2847
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human MCHR1 (NP_005288.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99705
Protein Weight	46kDa
Shipping	Ice bag