


Anti-Mcm4 Monoclonal Antibody For Wb, Elisa - P33991

Item Number: CM0010611



Introduction

Rabbit monoclonal antibody targeting human MCM4 for Western blot and ELISA. Cross-reacts with mouse and rat MCM4. Ideal for DNA replication and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	MCM4 Rabbit mAb
Remarks/Alias	NKCD; CDC21; CDC54; IMD54; NKGCD; hCdc21; P1-CDC21; MCM4
Species	Human
GeneID (human)	4173
GeneID	4173
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 764-863 of human MCM4 (P33991)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33991
Protein Weight	97kDa
Shipping	Ice bag