


Anti-Me1 Polyclonal Antibody (Wb, Elisa) - P48163

Item Number: CM0027785



Introduction

Rabbit polyclonal antibody against human ME1 (NADP-dependent malic enzyme). Validated for Western blot and ELISA, cross-reacts with mouse and rat. Detects a band at ~64 kDa. Suitable for metabolic research.

[Learn More](#)

Parameter	Details
Product Name	ME1 Rabbit pAb
Remarks/Alias	MES; HUMNDME; ME1
Species	Human
Gene ID (Human)	4199
Gene ID	4199
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ME1 (NP_002386.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48163
Protein Weight	64kDa
Shipping	Ice bag