

Anti-Mecp2 Polyclonal Antibody For Wb, Elisa - P51608

Item Number: CM0027033



Introduction

Anti-MECP2 Rabbit Polyclonal Antibody (P51608) validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Suitable for chromatin biology, DNA methylation, and neurodevelopmental disorder studies.

[Learn More](#)

Parameter	Value
Product Name	MECP2 Rabbit pAb
Remarks/Alias	RS; RTS; RTT; PPMX; MRX16; MRX79; MRXSL; AUTSX3; MRXS13; MECP2
Species	Human
GeneID (human)	4204
GeneID	4204
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 387-486 of human MECP2 (NP_004983.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51608
Protein Weight	52kDa
Shipping	Ice bag