

Anti-Med4 Polyclonal Antibody For Wb, Elisa - Q9Npj6

Item Number: CM0026071



Introduction

Rabbit polyclonal antibody against human MED4, validated for Western blot (WB) and ELISA. MED4 is a subunit of the Mediator complex, bridging regulatory proteins and RNA polymerase II for transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	MED4 Rabbit pAb
Remarks/Alias	ARC36; VDRIP; DRIP36; TRAP36; HSPC126; MED4
Species	Human
GeneID (Human)	29079
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human MED4 (NP_054885.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NPJ6
Protein Weight	30kDa
Shipping	Ice bag