


Anti-Mef2D Polyclonal Antibody For Wb, Ihc, Elisa - Q14814

Item Number: CM0019142



Introduction

Rabbit polyclonal antibody targeting human MEF2D (56 kDa), validated for WB, IHC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying transcriptional regulation in muscle and neuronal development.

[Learn More](#)

Parameter	Value
Product Name	MEF2D Rabbit pAb
Remarks/Alias	MEF2D
Species	Human
GeneID (Human)	4209
GeneID	4209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 431-521 of human MEF2D (NP_005911.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14814
Protein Weight	56kDa
Shipping	Ice bag