


Anti-Mek7 Rabbit Polyclonal Antibody (Ko Validated) For Wb, If/Icc, Elisa - O14733

Item Number: CM0020674



Introduction

KO validated rabbit polyclonal antibody targeting MEK7 (MKK7, JNKK2). Suitable for WB, IF/ICC, and ELISA in human, mouse, rat samples. Ideal for stress-activated JNK MAP kinase pathway and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MEK7 Rabbit pAb
Remarks/Alias	MEK; MKK7; JNKK2; MEK 7; MAPKK7; PRKMK7; SAPKK4; SAPKK-4; K7
Species	Human
GeneID (human)	5609
GeneID	5609
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human MEK7 (NP_660186.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14733
Protein Weight	47kDa
Shipping	Ice bag