

Anti-Mepce Rabbit Polyclonal Antibody For If/Icc And Elisa - Q7L2J0

Item Number: CM0029771



Introduction

Rabbit polyclonal antibody against human MEPCE (7SK snRNA methylphosphate capping enzyme), suitable for IF/ICC and ELISA. Ideal for studying 7SK RNP complex and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	MEPCE Rabbit pAb
Remarks/Alias	BCDIN3; MEPCE
Species	Human
Gene ID (Human)	56257
Gene ID	56257
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-689 of human MEPCE (NP_062552.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q7L2J0
Protein Weight	74kDa
Shipping	Ice bag