


Anti-Mettl13 Polyclonal Antibody For Wb, Elisa - Q8N6R0

Item Number: CM0030870



Introduction

Rabbit polyclonal antibody against human METTL13, validated for WB and ELISA. Detects METTL13 in human, mouse, and rat samples. Suitable for studying eEF1A methylation and translational regulation. Protein weight ~79kDa.

[Learn More](#)

Parameter	Value
Product Name	METTL13 Rabbit pAb
Remarks/Alias	feat; DFNM1; CGI-01; DFNB26; DFNB26M; KIAA0859; EEF1AKNMT; 5630401D24Rik; METTL13
Species	Human
GeneID (Human)	51603
GeneID	51603
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 540-699 of human METTL13 (NP_057019.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6R0
Protein Weight	79kDa
Shipping	Ice bag