

Anti-Mettl3 Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q86U44

Item Number: CM0021077



Introduction

High-quality rabbit polyclonal antibody against human METTL3 (M6A; IME4; Spo8; MT-A70; hMETTL3). Validated for WB, IHC-P, IP, ELISA. Cross-reacts with human, mouse, rat. For epitranscriptomics, m6A methylation, cancer research.

[Learn More](#)

Parameter	Value
Product Name	METTL3 Rabbit pAb
Remarks/Alias	M6A; IME4; Spo8; MT-A70; hMETTL3; METTL3
Species	Human
GeneID (Human)	56339
GeneID	56339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human METTL3 (NP_062826.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86U44
Protein Weight	64kDa
Shipping	Ice bag