


Anti-Mif Polyclonal Antibody For Wb, If/Icc, Elisa - P14174

Item Number: CM0013918



Introduction

Rabbit polyclonal antibody against human MIF (macrophage migration inhibitory factor), validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for innate immunity and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	MIF Rabbit pAb
Remarks/Alias	GIF; GLIF; MMIF; MIF
Species	Human
Gene ID (Human)	4282
Gene ID	4282
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human MIF (NP_002406.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14174
Protein Weight	12kDa
Shipping	Ice bag