

Anti-Mld Monoclonal Antibody For Wb, Elisa - 015121

Item Number: CM0003347

Introduction

Rabbit monoclonal antibody targeting human sphingolipid delta(4)-desaturase DES1 (MLD/DEGS1) for WB and ELISA. Ideal for lipid metabolism and leukodystrophy research. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	MLD Rabbit mAb
Remarks/Alias	MLD; DEGS; DES1; Des-1; FADS7; HLD18; MIG15; DEGS-1
Species	Human
GeneID (Human)	8560
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human MLD (NP_003667.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15121
Protein Weight	38kDa
Shipping	Ice bag