


Anti-Mlh1 Monoclonal Antibody (Ko Validated) For Wb, Elisa - P40692

Item Number: CM0007159



Introduction

KO validated mouse monoclonal antibody against human DNA mismatch repair protein Mlh1 (MLH1), suitable for Western blot and ELISA. Targets a key MMR component implicated in colorectal cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MLH1 Mouse mAb
Remarks/Alias	FCC2; COCA2; HNPCC; MLH-1; hMLH1; HNPCC2; LYNCH2; MMRC51; H1
Species	Mouse
GeneID	4292
Immunogen	Recombinant protein of human MLH1
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P40692
Protein Weight	85kDa
Shipping	Ice bag