


Anti-MLst8 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bvc4

Item Number: CM0014603



Introduction

Rabbit polyclonal antibody targeting MLST8 (GbetaL/mLST8), a subunit of mTORC1 and mTORC2; validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Ideal for nutrient signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MLST8 Rabbit pAb
Remarks/Alias	GBL; LST8; POP3; WAT1; GbetaL; MLST8
Species	Human
Gene ID	64223
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human MLST8 (NP_001186103.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BVC4
Protein Weight	36kDa
Shipping	Ice bag